

Summary Results: 2026-10-07_P9REAL_click_grid_run1

Subject: P9REAL · Test Date: Oct 7, 2026 5:09 PM · Evaluator: —

Test Configuration

Configuration Name: Standard

Setting	Value
Configuration name	Standard
Task	Grid Click
Input	Eye gaze (dwell), GP3HD, 150 Hz
Trials (planned)	18
Selection	Dwell 0.8 s, refractory 0.5 s
Target size	Medium (5°), 103 px radius
Layout	3×3 grid, gap Standard
Maximum time per trial	8 s
Pause between trials	0.8 s
Theme	Forest
Gaze cursor shown	Yes
Feedback	Sound on, sparkle on
Gaze smoothing	Smoothing α 0.22, jitter tolerance 40 px
Display	1920×1080 @ 100% (standard)
Viewing distance	650 mm
Calibration	5 points, mean error 73.7 px (1.78°), measured
Canvas	1920×965 physical px

Notes

One cell of a visible board lights up; the child selects it by looking at it.

Summary of Results

	% (N)	Trial Time (s)	Reaction Time (s)	Entries
Error-free Target Selections	44% (8/18)	1.38	0.57	1.0
All Targets Selected	100% (18/18)	1.86	0.53	1.8
Targets Not Selected	0% (0/18)	—	—	—
All Trials	100% (18/18)	1.86	0.53	1.8

Target area = drawn target + 40 px tolerance ring. Reaction Time = onset to the first gaze entry; about 0 if the gaze already rested on the new target's place.

Target Map



Eye Metrics

Metric	Value
Fixations	124 (6.9 per trial)
Mean fixation duration	0.22 s (median 0.19 s)
Saccades	568
Mean saccade amplitude	1.9°
Mean peak saccade velocity	136 deg/s (max 1623)
Scan-path length per trial	59.7°
Mean pupil diameter	4.10 mm
Mean pupil change from baseline	+0.01 mm (+0.4 %)
Valid gaze during trials	96 %
Calibration error	1.8° (74 px), measured

Trial-by-Trial Results

Trial	Size (deg)	Distance (deg)	Outcome	Trial Time (s)	Reaction Time (s)	Entries	Fixations	Mean fix. dur. (s)	Saccades	Mean peak vel. (deg/s)	Pupil (mm)	Pupil change (mm)
1	5.0	—	Hit	1.41	0.60	1	8	0.16	32	212	4.18	+0.14
2	5.0	13.1	Hit	5.44	0.31	4	20	0.19	79	160	3.88	-0.31
3	5.0	11.8	Hit	1.40	0.32	2	5	0.25	21	104	4.13	+0.50
4	5.0	16.6	Hit	1.12	0.32	1	4	0.23	16	140	3.92	-0.30
5	5.0	25.9	Hit	1.87	1.07	1	6	0.26	29	120	4.13	-0.01
6	5.0	23.2	Hit	2.07	0.63	2	9	0.21	49	120	4.37	+0.19
7	5.0	25.9	Hit	1.44	0.63	1	6	0.23	26	124	3.77	—
8	5.0	5.9	Hit	2.75	0.97	2	8	0.31	34	134	4.15	+0.28
9	5.0	11.7	Hit	1.18	0.01	2	3	0.24	20	109	4.02	-0.08
10	5.0	13.1	Hit	2.41	1.25	2	6	0.34	35	110	3.98	+0.06

Trial	Size (deg)	Distance (deg)	Outcome	Trial Time (s)	Reaction Time (s)	Entries	Fixations	Mean fix. dur. (s)	Saccades	Mean peak vel. (deg/s)	Pupil (mm)	Pupil change (mm)
11	5.0	23.9	Hit	1.37	0.26	2	5	0.22	20	89	3.92	+0.14
12	5.0	5.9	Hit	1.18	0.37	1	4	0.22	21	177	3.98	-0.09
13	5.0	11.8	Hit	3.12	0.53	5	16	0.17	64	138	4.44	+0.25
14	5.0	16.6	Hit	1.18	0.38	1	6	0.14	19	179	3.80	-0.62
15	5.0	11.8	Hit	1.32	0.52	1	5	0.23	24	107	4.25	+0.12
16	5.0	11.7	Hit	1.52	0.72	1	5	0.25	33	148	4.38	-0.00
17	5.0	5.9	Hit	1.52	0.46	2	5	0.25	22	121	4.27	-0.14
18	5.0	11.7	Hit	1.20	0.17	2	3	0.33	24	107	4.16	+0.01

Definitions

- Trial Time: from the target appearing to its selection (Hit) or to the end of the trial (Not selected).
- Reaction Time: from the target appearing to the first time the gaze entered the target area; about 0 if the gaze already rested on the new target's place.
- Entries: how many separate times the gaze entered the target area (leaving for less than 0.12 s does not count as a new entry).
- Error-free Target Selection: a Hit selected on the first gaze entry, with no earlier attempt.
- Fixations follow the Gazepoint definition and count for the trial in which they start.
- Saccades are found by velocity threshold (50 deg/s) on the device-rate gaze; peak velocity is a smoothed value, comparable only between children measured on the same device and sample rate.
- Pupil change: the mean diameter during the trial minus the mean of the 0.3 s before the target appeared; blinks (0.1 s either side of any invalid sample) are excluded.
- Skipped trials are excluded from every percentage.